

MOVIPREP

THIS PREP REQUIRES A PRESCRIPTION FROM YOUR DOCTORS OFFICE

THE MOVIPREP CARTON CONTAINS 4 POUCHES AND A DISPOSABLE CONTAINER FOR MIXING. YOU MUST COMPLETE THE ENTIRE PREP TO ENSURE THE MOST EFFECTIVE CLEANSING.

STEP 1

BEGINNING AT: 2:00 PM

EMPTY 1 POUCH A AND 1 POUCH B INTO THE DISPOSABLE CONTAINER.

ADD LUKE WARM DRINKING WATER TO THE TOP LINE OF THE CONTAINER.
MIX TO DISSOLVE

STEP 2

THE MOVIPREP CONTAINER IS DIVIDED BY 4 MARKS.
EVERY 15 MINUTES, DRINK THE SOLUTION DOWN TO THE NEXT MARK
(APPROXIMATELY 8 OZ), UNTIL THE FULL LITER IS COMPLETE

DRINK 16 OZ OF THE CLEAR LIQUID OF YOUR CHOICE.

BEGINNING AT: 8:00 PM

REPEAT STEPS 1 AND 2.

**IF YOU HAVE ANY QUESTIONS, PLEASE CALL OUR OFFICE AT
(310) 539-2055 DIAL 0 AND ASK FOR:**

CATRINA WILMA ISELA SUSIE JODIE